

IRRIGATOR™

IrriGator™ Plus, Pro, and Premium Irrigation Tips deliver the same clinical performance as other top rated irrigating tips, but at a *substantially* lower cost.



PREMIUM

Pre-Bent Dual Port Irrigating Tips

- Same performance as **NaviTip** irrigating tips for **45% LESS!**
- Dual side ports deliver irrigants safely; which minimizes the potential for chemicals being expressed past the apex of the canal.
- 30-gauge cannula allows for easier navigation through canals
- Tips are bendable at the pre-bent point for easier access to canals
- Available in 22mm & 26mm lengths

\$49.99

Pack of 50



PRO

Straight 25mm Bendable Dual Port Irrigating Tips

- Same performance as **Endo-Eze** irrigating tips for **40% LESS!**
- Luer lock hub threads to lock tip in place for increased safety.
- Closed-end tip for side-port reducing risk of expressing chemicals past the apex.
- Provides ideal reach in the canal.
- Flexible, cannula with an anti-obturing end.
- Available in 27-gauge and 30-gauge sizes.
- Tips come individually wrapped and pre-sterilized.

\$59.99

Pack of 100



PLUS

Straight 25mm Bendable Single Port Irrigating Tips

- Same performance as **Max-i-Probe**, **ProRinse**, and **Vista-Probe** tips at **45% LESS!**
- Luer lock hub threads to lock tip in place for increased safety.
- Closed-end tip for side-port reducing risk of expressing chemicals past the apex.
- Provides ideal reach in the canal.
- Flexible, cannula with an anti-obturing end.
- Available in 27-gauge and 30-gauge sizes.
- Tips come individually wrapped and pre-sterilized.

\$49.99

Pack of 100



NaviTip and Endo-Eze are registered trademarks of Ultradent Products, Inc. Max-i-Probe and ProRinse are registered trademarks of Dentsply International, Inc. Vista-Probe is a registered trademark of Vista Apex Dental Products.